



LOMAS /
MONDEFIT

The most innovative anchorage method for tooth correction

- Controlled and reliable tooth movement
- Completion of tooth correction without gaps
- For the young and young at heart







The most innovative anchorage method for tooth correction

The MONDEFIT system was developed based on the LOMAS system. Due to the compatibility of both systems by using the same instruments, this combination covers many orthodontic indications with amazing results.



MONDEFIT / ORTHODONTICS

The MONDEFIT system can be used for a variety of applications in the maxilla because of the screw design and the abutments available. The molar anchorage, the distal extension of the dental arch as well as the mesialization of molars can be named as important indications.

LOMAS / ORTHODONTICS

The usage of Mini-Anchor-Screws has become a standard for orthodontic treatment. The self-drilling LOMAS screws offer many advantages such as high stability, immediate screw loadability as well as reliable and secure fixation. There are additional applications by combining the LOMAS system with the innovative MONDEFIT system.



Indications and Advantages

Possible Indications



- Mesialization / Distalization of molars
- Molar uprighting
- Maximum anchorage (Front anchorage / Molar anchorage)
- Temporary dentures
- Maxillary expansion
- Opening gaps in case of retention
- Symmetrical alignment of asymmetries
- Displaced teeth


Advantages



- Immediate loading is possible
- Minimally invasive insertion protocol
- Orthodontists can insert by themselves
- Easy to use
- Easy to combine with applications of the orthodontic accessories
- Reduced treatment time



MONDEFIT - Screws (Packaging 1/ea.) Material: Titanium



MONDEFIT Screw, Ø 1.5 mm self-drilling

Art. No.	Diameter (d) x Length
A33-M4A-1507	1.5 x 7 mm
A33-M4A-1509	1.5 x 9 mm

MONDEFIT Screw, Ø 2.0 mm self-drilling

Art. No.	Diameter (d) x Length
A33-M4A-2007	2.0 x 7 mm
A33-M4A-2009	2.0 x 9 mm
A33-M4A-2011	2.0 x 11 mm

MONDEFIT Screw, Ø 1.5 mm self-drilling, sterile

Art. No.	Diameter (d) x Length
ST-A33-M4A-1507	1.5 x 7 mm
ST-A33-M4A-1509	1.5 x 9 mm

MONDEFIT Screw, Ø 2.0 mm self-drilling, sterile

Art. No.	Diameter (d) x Length
ST-A33-M4A-2007	2.0 x 7 mm
ST-A33-M4A-2009	2.0 x 9 mm
ST-A33-M4A-2011	2.0 x 11 mm



MONDEFIT Emerg. Screw, Ø 2.3 mm self-drilling

Art. No.	Diameter (d) x Length
A33-M4A-2309	2.3 x 9 mm

MONDEFIT Emerg. Screw, Ø 2.3 mm self-drilling, sterile

Art. No.	Diameter (d) x Length
ST-A33-M4A-2309	2.3 x 9 mm



MONDEFIT - Abutments

Material: Implant steel



A33-ME0-0001 Abutment Standard, incl. Abutment screw



A33-ME2-0001 Abutment with 1 bracket, incl. Abutment screw



A33-ME2-0002 Abutment with 2 brackets, incl. Abutment screw



A33-ME1-0002 Abutment with wire 1.1 mm (Length: 120 mm), incl. Abutment screw

MONDEFIT - Impression

Material: Steel



A33-VXI-0001 Impression Cap Metal



A33-VXI-0002 Laboratory Analog



MONDEPLATES Material: Implant steel



A33-MR0-0001 MONDEPLATE connection plate 1.2 mm, long (17.50 mm), incl. fixation screws



A33-MP0-0001 MONDEPLATE connection plate 1.2 mm, short (12.50 mm), incl. fixation screws



A33-MPI-1211 MONDEPLATE connection plate, 1.2 mm, short, with 1.1 mm wire (Length: 120 mm), incl. fixation screws

A33-MRI-1211 MONDEPLATE connection plate, 1.2 mm, long, with 1.1 mm wire (Length: 120 mm), incl. fixation screws



A33-MB0-0001 Spare fixation screws (2/ea.)



A33-MB0-0002 Spare fixation screws (2/ea.)

MONDEFIT - Accessories



A33-MZ2-0007 Hook lock (1/ea.)



A33-MZ2-0008 Mobilizer for wires (1/ea.)



A33-MZ2-0010 Hook lock extended (1/ea.)



A33-MZ2-0011 Mobilizer for wires, extended (1/ea.)



855-606 Spring 240g (10/ea.)






855-180 Coil spring 150g (10/ea.)

58001182 Spring 500g (10/ea.)

855-179 NiTi Coil spring 6mm 200g (10/ea.)



LOMAS Screws (Packaging 1/ea.)



STANDARD Screw, \varnothing 1.5 mm self-drilling

HOOK Screw, \varnothing 1.5 mm self-drilling

Art. No.	Diameter (d) x Length
A33-11507	1.5 x 7 mm
A33-11509	1.5 x 9 mm
A33-11511	1.5 x 11 mm

Art. No.	Diameter (d) x Length
A33-21507	1.5 x 7 mm
A33-21509	1.5 x 9 mm
A33-21511	1.5 x 11 mm



STANDARD Screw, \varnothing 2.0 mm self-drilling

HOOK Screw, \varnothing 2.0 mm self-drilling

Art. No.	Diameter (d) x Length
A33-12207	2.0 x 7 mm
A33-12209	2.0 x 9 mm
A33-12211	2.0 x 11 mm

Art. No.	Diameter (d) x Length
A33-22007	2.0 x 7 mm
A33-22009	2.0 x 9 mm
A33-22011	2.0 x 11 mm

STANDARD Emerg. Screw, \varnothing 2.3 mm self-drilling


HOOK Emerg. Screw, \varnothing 2.3 mm self-drilling

Art. No.	Diameter (d) x Length
A33-12309	2.3 x 9 mm

Art. No.	Diameter (d) x Length
A33-22109	2.3 x 9 mm



LOMAS Screws (Packaging 1/ea.)



QUAD Screw H, Ø 1.5 mm
q = 0.018 x 0.025", self-drilling

Art. No.	Diameter (d) x Length
A33-31807	1.5 x 7 mm
A33-31809	1.5 x 9 mm
A33-31811	1.5 x 11 mm




QUAD Screw H, Ø 1.5 mm
q = 0.022 x 0.028", self-drilling

Art. No.	Diameter (d) x Length
A33-32207	1.5 x 7 mm
A33-32209	1.5 x 9 mm
A33-32211	1.5 x 11 mm



QUAD Screw H, Ø 2.0 mm
q = 0.018 x 0.025", self-drilling

Art. No.	Diameter (d) x Length
A33-41807	2.0 x 7 mm
A33-41809	2.0 x 9 mm
A33-41811	2.0 x 11 mm



QUAD Screw H, Ø 2.0 mm
q = 0.022 x 0.028", self-drilling

Art. No.	Diameter (d) x Length
A33-42207	2.0 x 7 mm
A33-42209	2.0 x 9 mm
A33-42211	2.0 x 11 mm

QUAD Emerg. Screw H, Ø 2.3 mm
q = 0.018 x 0.025", self-drilling



Art. No.	Diameter (d) x Length
A33-41909	2.3 x 9 mm

QUAD Emerg. Screw H, Ø 2.3 mm
q = 0.022 x 0.028", self-drilling

Art. No.	Diameter (d) x Length
A33-42309	2.3 x 9 mm




LOMAS Screws (Packaging 1/ea.)



suitable for arch wires

V=vertical



QUAD Screw V, Ø 2.0 mm
q = 0.018 x 0.025", self-drilling

QUAD Screw V, Ø 2.0 mm
q = 0.022 x 0.028", self-drilling

Art. No.	Diameter (d) x Length
A33-51807	2.0 x 7 mm
A33-51809	2.0 x 9 mm
A33-51811	2.0 x 11 mm

Art. No.	Diameter (d) x Length
A33-52207	2.0 x 7 mm
A33-52209	2.0 x 9 mm
A33-52211	2.0 x 11 mm

QUAD Emerg. Screw V, Ø 2.3 mm
q = 0.018 x 0.025", self-drilling


QUAD Emerg. Screw V, Ø 2.3 mm
q = 0.022 x 0.028", self-drilling

Art. No.	Diameter (d) x Length
A33-51909	2.3 x 9 mm

Art. No.	Diameter (d) x Length
A33-52309	2.3 x 9 mm



LOMAS VEGAS Screws (Packaging 1/ea.)



VEGAS Screw H, Ø 2.0 mm
q = 0.022 x 0.028", self-drilling, sterile

Art. No.	Diameter (d) x Length x Thread
ST-A33-53209	2.0 x 9 mm x 7 mm
ST-A33-53211	2.0 x 11 mm x 7.5 mm
ST-A33-53213	2.0 x 13 mm x 8 mm
ST-A33-53215	2.0 x 15 mm x 8.5 mm
ST-A33-53217	2.0 x 17 mm x 9 mm
ST-A33-53220	2.0 x 20 mm x 9.5 mm
ST-A33-53225	2.0 x 25 mm x 10 mm
ST-A33-53230	2.0 x 30 mm x 10.5 mm
ST-A33-53235	2.0 x 35 mm x 11 mm
ST-A33-53240	2.0 x 40 mm x 11.5 mm



Instruments

ML For MONDEFIT, LOMAS and LOMAS VEGAS:



A33-18232 Blade for angled handpiece, short, 23 mm, DENTAL

A33-18237 Blade for angled handpiece, long, 33 mm, DENTAL



A12-12625 Drill 1.0 x 30 mm, WL 15 mm, DENTAL
for 1.5 mm screws



A10-67517 Drill 1.5 x 28 mm, WL 15 mm, DENTAL
for 2.0 mm screws



A10-61000 Screwdriver handle, rotatable, 9.5 cm



A33-18230 Blade, long for A10-61000



A33-SS4-1100 Screwdriver, fixed



Instruments

M For MONDEFIT only:



A84-SK0-2000 MCD blade, 20 mm, DENTAL, self-retaining, for angled handpiece



A33-SS0-1000 Handle with blade

For fixation of:

- Fixation screws
- Abutments

For activating:

- Mobilizer



A99-003-04 Screwdriver handle, short, 6 cm (for laboratory use only)



A99-003-05  MCD blade, self-retaining (for laboratory use only)

LV For LOMAS VEGAS only:



A85-VA2-1500 Drill 1.5 x 60 mm, WL 15 mm, DENTAL for LOMAS VEGAS-Screws



Container

ML For MONDEFIT and LOMAS:



A73-31960 Sterilizing tray for implants and instruments, empty, incl. lid and length clips



A75-TC3-0677-7	Length clip 7 mm
A75-TC3-0677-9	Length clip 9 mm
A75-TC3-0677-11	Length clip 11 mm



MONDEFIT - Models

Ordering information



A95-13012
Model Mesialization / Distalization



A95-13007
Model Front anchorage



A95-13010
Model Maxillary expansion



A95-13002
Model Molar anchorage



Clinical case: Maximum anchorage



Maximum anchorage of the molars



Bilateral mesial space closure achieved by indirect anchorage

Clinical case: Distalization I



Initial situation



Bilateral tooth movement distally



Treatment outcome



Clinical case: Distalization 2



The distalization based on skeletal anchorage is evenly spreaded since the bicuspids (dentes premolares) move distal with the molars through the traction of the interdental ligaments. As a result, there are no large gaps.




Clinical case: Distalization 3






Clinical case: Mesialization



Mesialization of the upper molars through unilateral protraction palatinal.




Top 4 Areas of Application



- 1 Infrazygomatic
- 2 Interdental
- 3 Mandibular
- 4 Palatinal

Instruction for insertion of the screw





Clinical case: Molar uprighting and distalization



Insert LOMAS STANDARD screw as skeletal anchorage incident to the tooth using 2 NiTi coil springs



Position of the screw parallel to the tooth axis




Advanced uprighting and distalization



Completed uprighting and distalization



Placing the second dental implant



Placing the cap on the dental implants. After the screw has been removed, the orthodontic treatment is successfully completed.

Clinical case: Protraction of molars



Insertion of the LOMAS HOOK screw between the first and the second molar (interdental)



Transmission through rubber band




Transmission through NiTi coil spring

Progress of protraction →


Clinical case: Skeletal Gummy Smile Correction



Patient with Gummy Smile, overbite and tooth gap (diastema)



Positioning and function of the LOMAS QUAD screws



Transmission through rubber band as a tool for the correction of the overbite and for gap closure



Transmission through rubber band as a tool for the correction of the Gummy Smile




Transmission through rubber band as a tool for gap closure



Finished and activated appliance



Advanced correction front view



Advanced correction profile






Completed orthodontic treatment outcome



Features

LOMAS VEGAS Screw H, Ø 2.0 mm
 q = 0.022 x 0.028", self-drilling, sterile

Diameter (d) x Length x Thread
2.0 x 9 mm x 7 mm
2.0 x 11 mm x 7.5 mm
2.0 x 13 mm x 8 mm
2.0 x 15 mm x 8.5 mm
2.0 x 17 mm x 9 mm
2.0 x 20 mm x 9.5 mm
2.0 x 25 mm x 10 mm
2.0 x 30 mm x 10.5 mm
2.0 x 35 mm x 11 mm
2.0 x 40 mm x 11.5 mm



suitable for arch wires

H=horizontal


- Extended application areas: chin-menton, mandibular, (infra)zygomatic, interdental, palatal
- Alternative for existing orthodontic mini-plate-systems
- Fast healing around the long thread-free shank of the screw
- Opportunity of a stable anchorage in the mandible
- Sterile packaging




Developed in cooperation with: Dr. Villegas • specialist for Orthodontics and OMF surgeon • Colombia




Areas of Application



Zygomatic



Infrazygomatic



Interdental



Mandibular



Chin-menton



Clinical case: Class III patient



Original image



Class III, 11 years old patient




Original image



X-Ray view



Drilling in the maxilla



Placement of 2 LOMAS VEGAS screws (length 35 mm) in the zygomatic bone through 1 cm incision



Drilling in the mandible



Flapless placement of 2 LOMAS VEGAS screws (length 30 mm) in the mandibular symphysis



Placed screws for the orthopedic treatment



Soft tissue healing around the screws in the mandible




Class III elastics for intraoral class III orthopedic treatment



Attachment of an auxiliary wire with a hook in the rectangular slot to facilitate the placement of the elastics on the screws in the zygoma




Clinical case: Temporary dentures



Patient with missing upper permanent lateral incisors and canines



Tomography cut showing bone availability and root resorption of the upper right deciduous canine



LOMAS VEGAS screw (length 20 mm) was placed flapless



LOMAS VEGAS in place



Composite was added to the head to serve as an abutment



Temporary prosthesis in place



Tomographic image of the screw




Notice the great soft tissue healing around the LOMAS VEGAS with the new temporary prosthesis after 10 months



Clinical case: Comparison LOMAS VEGAS vs. 3D-VESTIPLATE



Pre-surgical occlusion and profile of a class III skeletal patient



Surgery first approach was performed placing a 3D-VESTIPLATE with 2 mini screws on the right side ...




... leaving increased overbite ...



... and a long LOMAS VEGAS screw (length 40 mm) on the left zygoma to be used as skeletal anchorage during post-surgical orthodontic treatment




Placing 2 LOMAS VEGAS screws (length 20 mm) in the mandibular buccal bilaterally in order to have anchorage to distalize the lower dental arch




NiTi coil springs used to distalize (right side)




NiTi coil springs used to distalize (left side)



Soft tissue healing around the screws and 3D-VESTIPLATE after 5 months of treatment: Maxilla right



Soft tissue healing around the screws and 3D-VESTIPLATE after 5 months of treatment: Maxilla left



Soft tissue healing around the screws and 3D-VESTIPLATE after 5 months of treatment: Mandible right



Soft tissue healing around the screws and 3D-VESTIPLATE after 5 months of treatment: Mandible left



Final occlusion 13 months after surgery



Instruction: Handling of sterile screws



Sterile version of the MONDEFIT, LOMAS and LOMAS VEGAS Screws



Pushing in the cardboard packaging above the perforation




Opening of the cardboard packaging and removing of the sterile double blister packaging



Sterile double blister packaging



Opening of the first blister and removing of the second sterile blister




Second sterile blister




Opening of the second sterile blister



Attention: Before folding back the blister to expose the screw head, make sure the screw head is not pulled out since especially in case of using the short screws (7-9 mm) it could happen that they fall out of the tube. Therefore: Do not put your thumb over the screw head.



Folding back the blister to expose the screw head

**Instruction: Handling of sterile screws**

Receiving the screw with the screwdriver / the blade
(horizontal and with little pressure)



The screw sticks in the blade.
Control: Only 1-2 tenths of the screw
head platform should be seen.

Applicable for these products:

- MONDEFIT
- LOMAS
- LOMAS VEGAS



Orthodontics by AFM

We reserve the right to make catalog and design changes arising from further developments and modifications. Pictures, product descriptions and texts are the property of Arno Fritz GmbH. Further use and reproduction by Third Parties requires written consent from Arno Fritz GmbH. All rights reserved! Note: Responsibility for the selection of the patient, for adequate training as well as the decision as to the choice and use of the implants, a post-operative removal or retention of implants is the sole decision of the surgeon. All brands and logos referred to in this brochure, registered by Third Parties, are subject without restriction to the legal provisions currently in force and to the rights of ownership of the respective registered owner.

© 2016 - Arno Fritz GmbH



Arno Fritz GmbH • Am Gewerbering 7 • 78570 Mühlheim a. d. Donau/Germany
Phone +49 7463 99307 60 • Fax +49 7463 99307 66 • info@arno-fritz.com

