Scan for more information

Check out our new learning platform Reliance U! www.RelianceU.tech

Reliance Orthodontic Products, Inc.

1540 W. Thorndale Avenue, Itasca, IL 60143 USA

1-800-323-4348 | Made in USA | Rx Only
RelianceOrthodontics.com

EMERGO EUROPE

Prinssgracht 20, 2514 AP The Hague, The Netherlands

AUSTRALIAN SPONSOR

EMERGO AUSTRALIA

Level 20, Tower II, Darling Park, 201 Sussex St, Sydney, NSW 2000, Australia

1ST & FINAL®

V1 NOV.2021

RelianceOrthodontics.com

01-EN ENGLISH TRANSLATION

INDICATIONS FOR USE: Orthodontic Prophy Pumice for use in cleaning the tooth surface prior to acid etching.

PRODUCT INFORMATION:

- Pre-packaged and pre-mixed orthodontic prophy pumice.
- No flavoring, fluorides, oils or chemicals which could leave a film that may inhibit or weaken the bond.
- 1st & Final® prophy pumice can be used prior to bonding to achieve a thorough prophylaxis; or after appliance removal to achieve a clean fresh surface.

INSTRUCTIONS FOR USE

1. Remove foil barrier from 1st & Final® prophy cup.
 2. Using a prophy angle (Reliance Prep Step™ Disposable Prophy Angle), fill the prophy angle cup with paste and establish a firm finger rest.
 3. Bring the prophy cup to close proximity to the tooth surface before operating the prophy angle.
 4. Replenish the supply of paste in the prophy cup, as needed. Do not overfill.
 5. After each patient, dispose of prophy angle and unused prophy cup of 1st & Final®. Do Not Re-use.
- ② Single Use Only. Each individual unit of 1st & Final® contains sufficient paste for one full prophylaxis procedure.

⚠ Excessive prophylaxis paste should be expectorated, not swallowed. Will cause eye irritation. May cause skin irritation. Prolong exposure may cause respiratory system irritation. Ingestion of large amounts may cause nausea, vomiting, headaches, or dizziness. Refer to SDS for additional information.

WARRANTY: Reliance Orthodontic Products, Inc. recognizes its responsibility to replace products if proven to be defective. Reliance Orthodontic Products, Inc. does not accept liability for any damages or loss, direct or consequential, stemming from the use of or inability to use the products as described. Before using, it is the responsibility of the user to determine the suitability of the product for its intended use. The user assumes all risk and liability in connection therewith.

🌡️ **TEMPERATURE LIMITATION:** Store device at 15-30 Degrees C).

RX ONLY: Federal law restricts this device to sale by or on the order of a dental professional.

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GARANTIE: Reliance Orthodontic Products, Inc. își recunoaște responsabilitatea de a înlocui produsele