



GC Ortho Connect™ Flow

Light Cure Orthodontic Adhesive
One Step System



Since 1921
Towards Century of Health

Experts in Orthodontics
GC Orthodontics Europe GmbH
www.gcorthodontics.eu




GC Ortho **Connect™** Flow 90-3100-0000 (2 syringes of 2,7 g each)

Light Cure Orthodontic Adhesive for lingual retainers, stops and aligners attachments.

Easy to handle

Like GC Ortho Connect, GC Ortho Connect flow contain the primer in the paste so you save one step in the bonding process. Easy to fixing lingual retainer wire.



Clean and dry




Etch





Set the retainer



Apply Ortho Connect Flow



Light Cure




Improved flexibility and toughness


Thanks to a new monomer with long chain flexible, the flexibility of the material is enhanced.

High Flexibility: High retention force to retainer wire

High Toughness: The resin is hard to break



GC Ortho Connect Flow
GC Ortho Connect
Competitor 1
Competitor 2



Strong adhesion to the retainer wire or archwire

Ortho Connect Flow is self-adhesive and bonds well to both enamel and retainer wire. It contains MDP like metal primer which ensures a high bond strength to the wire.

GC Orthodontics Europe GmbH
Harkortstrasse 2
D-58339 Breckerfeld
Germany

Direction commerciale
8, rue Benjamin Franklin
94370 Sucy en Brie
Tél. : 01 77 45 65 61
Fax : 01 45 76 31 53
info.gco.france@gc.dental
www.gcorthodontics.eu

